



ISO 9001:2015 Certified

SAFETY DATA SHEET

according to Federal Register / Vol. 77, No. 58 / Monday, March 26, 2012 / Rules and Regulations

Preparation Date: 01/01/2019

Revision Date: N/A

Revision Number: N/A

01 Identification of the substance/mixture and of the company/undertaking

- Trade name : Fmoc-Glu(OtBu)-OH · H₂O
- Catalog : AA31283
- CAS Number: 204251-24-1
- Relevant identified uses of the substance or mixture and uses advised against
- Application of the substance / the preparation
Laboratory chemicals
- Details of the supplier of the safety data sheet
- Manufacturer/Supplier:
Dawn Scientific Inc
121 Liberty Street, Metuchen, NJ, 08840
Tel : 732-902-6300
Fax : 973-802-1005
sales@dawnscientific.com
www.dawnscientific.com
- Further information obtainable from:
Department: Marketing
- Emergency telephone number:
+41 44 251 51 51 (Tox Info Suisse)

*

02 Hazards identification

- Classification of the substance or mixture
- Classification according to Regulation (EC) No 1272/2008
Void
- Label elements
- Labelling according to Regulation (EC) No 1272/2008
Void
- Signal word
Void
- Hazard statements
Void
- Other hazards
- Results of PBT and vPvB assessment
- PBT:
Not applicable.
- vPvB:
Not applicable.

03 Composition/information on ingredients

- Chemical characterization: Substances
- | CAS No. | Description |
|----------------|--------------------------------------|
| 204251-24-1 | Fmoc-Glu(OtBu)-OH · H ₂ O |
- Identification number(s)
 - EINECS Number:
-
 - RTECS Number:
-

04 First aid measures

- General information:
No special measures required.
- After inhalation:
Supply fresh air; consult doctor in case of complaints.
- After skin contact:
Immediately wash with water and soap and rinse thoroughly.
- After eye contact:
Rinse opened eye for several minutes under running water.
- After swallowing:
Call a doctor immediately.
- Information for doctor:
- Most important symptoms and effects, both acute and delayed
No further relevant information available.
- Indication of any immediate medical attention and special treatment needed
No further relevant information available.

05 Firefighting measures

- Suitable extinguishing agents:
CO₂, powder or water spray. Fight larger fires with water spray or alcohol resistant foam.
- Special hazards arising from the substance or mixture
No further relevant information available.
- Advice for firefighters
- Protective equipment:
No special measures required.

06 Accidental release measures

- Personal precautions, protective equipment and emergency procedures
Not required.
- Environmental precautions:
No special measures required.
- Methods and material for containment and cleaning up:
Pick up mechanically.

07 Handling and storage

- Handling:
- Precautions for safe handling
Prevent formation of dust.
- Information about fire - and explosion protection:
No special measures required.
- Storage:
- Requirements to be met by storerooms and receptacles:
No special requirements.
- Information about storage in one common storage facility:
Store away from foodstuffs.
- Further information about storage conditions:
Store in dry conditions.
- Recommended storage temperature:
+5 +/- 3 °C
- Specific end use(s)
No further relevant information available.

08 Exposure controls/personal protection

- Ingredients with limit values that require monitoring at the workplace:
Not required.
- Additional information:
The lists valid during the making were used as basis.
- Personal protective equipment:
- General protective and hygienic measures:
The usual precautionary measures are to be adhered to when handling chemicals.
- Respiratory protection:
In case of brief exposure or low pollution use respiratory filter device. In case of intensive or longer exposure use self-contained respiratory protective device.
Use suitable respiratory protective device in case of insufficient ventilation.
- Protection of hands:
The glove material has to be impermeable and resistant to the product/ the substance/ the preparation.
Due to missing tests no recommendation to the glove material can be given for the product/ the preparation/ the chemical mixture.
Selection of the glove material on consideration of the penetration times, rates of diffusion and the degradation
Protective gloves
- Material of gloves
The selection of the suitable gloves does not only depend on the material, but also on further marks of quality and varies from manufacturer to manufacturer.
- Penetration time of glove material
The exact break through time has to be found out by the manufacturer of the protective gloves and has to be observed.
- Eye protection:
Safety glasses

09 Physical and chemical properties

| | |
|--|--|
| General Information | |
| Appearance: | |
| Form: | Powder |
| Colour: | White |
| Odour: | Nearly odourless |
| Odour threshold: | Not determined. |
| Molecular Weight | 443,5000 g/mol |
| Change in condition | |
| Melting point/Melting range: | Undetermined. |
| Boiling point/Boiling range: | Undetermined. |
| Flash point: | Not applicable. |
| Flammability (solid, gaseous): | Not applicable. |
| Ignition temperature: | |
| Decomposition temperature: | Not determined. |
| Self-igniting: | Not determined. |
| Danger of explosion: | Not determined. |
| Explosion limits: | |
| Lower: | Not determined. |
| Upper: | Not determined. |
| Vapour pressure: | Not determined. |
| Density: | Not determined. |
| Solubility in / Miscibility with water: | Not determined. |
| Viscosity: | |
| Dynamic: | Not determined. |
| Kinematic: | Not determined. |
| Other information | No further relevant information available. |

10 Stability and reactivity

- Thermal decomposition / conditions to be avoided:
No decomposition if used according to specifications.
- Possibility of hazardous reactions
No dangerous reactions known.
- Conditions to avoid
No further relevant information available.
- Incompatible materials:
No further relevant information available.
- Dangerous decomposition products:
No dangerous decomposition products known.

11 Toxicological information

- Acute toxicity:
No further relevant information available.
- Primary irritant effect:
- on the skin:
No further relevant information available.
- on the eye:
No further relevant information available.
- Sensitization:
No sensitizing effects known.
- Additional toxicological information:
The substance is not subject to classification according to the latest version of the EU lists.

12 Ecological information

- Aquatic toxicity:
No further relevant information available.
- Persistence and degradability
No further relevant information available.
- Behaviour in environmental systems:
- Bioaccumulative potential
No further relevant information available.
- Mobility in soil
No further relevant information available.
- Additional ecological information:
- General notes:
Not known to be hazardous to water.
- Results of PBT and vPvB assessment
- PBT:
Not applicable.
- vPvB:
Not applicable.
- Other adverse effects
No further relevant information available.

13 Disposal considerations

- Waste treatment methods
- Waste disposal key:
18 02 06
- Uncleaned packaging:
- Recommendation:
Disposal must be made according to official regulations.

14 Transport information

- | | |
|-------------------------|---|
| ADR/RID class: | - |
| IMDG Class: | - |
| ICAO/IATA Class: | - |
- Environmental hazards:
Not applicable.
 - Transport in bulk according to Annex II of MARPOL73/78 and the IBC Code
Not applicable.

15 Regulatory information

- Chemical safety assessment:
A Chemical Safety Assessment has not been carried out.

16 Other information

This information is based on our present knowledge. However, this shall not constitute a guarantee for any specific product features and shall not establish a legally valid contractual relationship.

- Department issuing MSDS:
Department: Marketing
- Abbreviations and acronyms:
ADR: Accord européen sur le transport des marchandises dangereuses par Route (European Agreement concerning the International Carriage of Dangerous Goods by Road)
RID: Règlement international concernant le transport des marchandises dangereuses par chemin de fer (Regulations Concerning the International Transport of Dangerous Goods by Rail)
IMDG: International Maritime Code for Dangerous Goods
IATA: International Air Transport Association
ICAO: International Civil Aviation Organisation
EINECS: European Inventory of Existing Commercial Chemical Substances
CAS: Chemical Abstracts Service (division of the American Chemical Society)
RTECS: Registry of Toxic Effects of Chemical Substances
- * Data compared to the previous version altered.