



SAFETY DATA SHEET

according to Federal Register / Vol. 77, No. 58 / Monday, March 26, 2012 / Rules and Regulations

Preparation Date: 01/01/2019

Revision Date: N/A

Revision Number: 1

SECTION 1 Identification**1.1. Product Identifier**

Trade Name or Designation: Phosphate Buffered Saline, 10X, pH 8.0 8% (w/v) NaCl / 0.2% (w/v) KCl

Product Number: B16053

1.2. Recommended Use and Restrictions on Use

General Laboratory Reagent

1.3. Details of the Supplier of the Safety Data Sheet**Dawn Scientific Inc**

121 Liberty Street, Metuchen, NJ, 08840

Tel : 732-902-6300 | Fax : 973-802-1005

sales@dawnscientific.com | www.dawnscientific.com

1.4. Emergency Telephone Number (24 hours)

CHEMTREC (USA) 800-424-9300

CHEMTREC (International) 1+ 703-527-3887

SECTION 2: Hazard(s) Identification**2.1. Classification of the Substance or Mixture**

For the full text of the Hazard and Precautionary Statements listed below, see Section 16.

Hazard Class	Category	Hazard Statements	Precautionary Statements:
Eye Damage / Irritation	Category 2A	H319	P264, P280, P305+P351+P338, P337+P313

2.2. GHS Label Elements

Pictograms:



Signal Word:

Warning

Hazard Statements:

Hazard Number	Hazard Statement
H319	Causes serious eye irritation.

Precautionary Statements:

Precautionary Number	Precautionary Statement
P264	Wash arms, hands and face thoroughly after handling.
P280	Wear protective gloves and eye protection.
P305+P351+P338	IF IN EYES: Rinse cautiously with water for several minutes. Remove contact lenses, if present and easy to do. Continue rinsing.
P337+P313	If eye irritation persists: Get medical attention.

2.4. Hazards not Otherwise Classified or Covered by GHS

Data not available.

SECTION 3: Composition / Information on Ingredients

3.1. Components of Substance or Mixture

Chemical Name	Formula	Molecular Weight	CAS Number	Weight%
Water	H ₂ O	18.01 g/mol	7732-18-5	90.77
Sodium Chloride	NaCl	58.44 g/mol	7647-14-5	7.48
Sodium Phosphate Dibasic	Na ₂ HPO ₄	141.95 g/mol	7558-79-4	1.35
Potassium Dihydrogen Phosphate	KH ₂ PO ₄	136.08 g/mol	7778-77-0	0.22
Potassium Chloride	KCl	74.55 g/mol	7447-40-7	0.19

SECTION 4: First-Aid Measures

4.1. General First Aid Information

Eye Contact: IF IN EYES: Rinse cautiously with water for several minutes. Remove contact lenses, if present and easy to do. Continue rinsing. May cause irritation, redness, pain, and tearing.

Inhalation: Not expected to require first aid. If necessary, remove to fresh air.

Skin Contact: May cause irritation, redness, itching and pain, especially on damaged or broken skin.

Ingestion: Dilute with water or milk. Call a physician if necessary.

4.2. Most Important Symptoms and Effects, Acute and Delayed

Causes serious eye irritation. May cause mild irritation to the eyes and skin. Wash areas of contact with water. Call a physician if irritation develops. Does not present any significant health hazards. EYE CONTACT: May cause irritation, redness, pain, and tearing. SKIN CONTACT: May cause irritation, redness, itching and pain, especially on damaged or broken skin.

4.3. Medical Attention or Special Treatment Needed

Not expected to require special treatment.

SECTION 5: Fire-Fighting Measures

5.1. Extinguishing Media

Use any means suitable for extinguishing surrounding fire.

5.2. Specific Hazards Arising from the Substance or Mixture

Not considered to be a fire or explosion hazard.

5.3. Special Protective Equipment for Firefighters

Use protective clothing and breathing equipment appropriate for the surrounding fire.

SECTION 6: Accidental Release Measures

6.1. Personal Precautions, Protective Equipment and Emergency Procedures

Wear protective gloves and eye protection.

6.2. Cleanup and Containment Methods and Materials

Absorb with suitable material (vermiculite, clay, etc.) and dispose of in accordance with local regulations. Check with local agencies for the proper disposal of phosphate containing solutions.

SECTION 7: Handling and Storage

7.1. Precautions for Safe Handling and Storage Conditions

As with all chemicals, wash hands thoroughly after handling. Avoid contact with eyes and skin. Protect from freezing and physical damage.

SECTION 8: Exposure Controls / Personal Protection

8.1 Control Parameters

Chemical Name	Limit Type	Country	Exposure Limit	Information Source
Data not available.				

8.2. Exposure Controls

Engineering Controls: No specific controls are needed. Normal room ventilation is adequate.

Respiratory Protection: Normal room ventilation is adequate.

Skin Protection: Wear protective gloves and eye protection. Chemical resistant gloves.

Eye Protection: Wear protective gloves and eye protection. Safety glasses or goggles.

8.3. Personal Protective Equipment

Wear protective gloves and eye protection. Normal room ventilation is adequate. Chemical resistant gloves. Safety glasses or goggles.

SECTION 9: Physical and Chemical Properties

9.1. Basic Physical and Chemical Properties

Appearance: Colorless liquid

Physical State: Liquid

Odor: Data not available.

Odor Threshold: Data not available.

pH: 8

Melting/Freezing Point: Data not available.

Initial Boiling Point/Range: Approximately 100.1°C - Approximately 100.1°C

Flash Point: Data not available.

Evaporation Rate: Data not available.

Flammability: Data not available.

Flammability/Explosive Limits: Data not available.

Vapor Pressure: Data not available.

Vapor Density: Data not available.

Relative Density: 1.1

Solubility: Miscible

Partition Coefficient: Data not available.

Auto-Ignition Temperature: Data not available.

Decomposition Temperature: Data not available.

Viscosity: Data not available.

Explosive Properties: Data not available.

Oxidizing Properties: Data not available.

SECTION 10: Stability and Reactivity

10.1. Reactivity and Chemical Stability

Stable under normal conditions of use and storage.

10.2. Possibility of Hazardous Reactions

Data not available.

10.3. Conditions to Avoid and Incompatible Materials

Strong oxidizing agents, Lithium, Bromine Trifluoride.

10.4. Hazardous Decomposition Products

Will not occur.

SECTION 11: Toxicological Information

11.1. Information on Toxicological Effects

Acute Toxicity - Oral Exposure:

Not applicable.

Acute Toxicity - Dermal Exposure:

Not applicable.

Acute Toxicity - Inhalation Exposure:

Not applicable.

Acute Toxicity - Other Information:

LD50, Oral, Rat: (Sodium Chloride) 3000 mg/kg, details of toxic effects not reported other than lethal dose value. LD50, Dermal, Rabbit: (Sodium Chloride) >10 gm/kg. LD50, Oral, Rat: (Sodium Phosphate Dibasic) 17 gm/kg; LD50, Dermal, Rabbit: (Potassium Phosphate Monobasic) >4640 mg/kg; details of toxic effects not reported other than lethal dose value. LD50, Oral, Rat: (Potassium Chloride) 2600 mg/kg; Irritation: eye, rabbit (500mg/24 hr mild).

Skin Corrosion and Irritation:

Not applicable.

Serious Eye Damage and Irritation:

Causes serious eye irritation. Wash arms, hands and face thoroughly after handling. Wear protective gloves and eye protection. IF IN EYES: Rinse cautiously with water for several minutes. Remove contact lenses, if present and easy to do. Continue rinsing. If eye irritation persists: Get medical attention.

Respiratory Sensitization:

Not applicable.

Skin Sensitization:

Not applicable.

Germ Cell Mutagenicity:

Not applicable.

Carcinogenicity:

Not applicable.

Reproductive Toxicity:

Not applicable.

Specific Target Organ Toxicity from Single Exposure:

Not applicable.

Specific Target Organ Toxicity from Repeated Exposure:

Not applicable.

Aspiration Hazard:

Not applicable.

Additional Toxicology Information:

Data not available.

SECTION 12: Ecological Information

12.1. Ecotoxicity

Not applicable.

12.2. Persistence and Degradability

Data not available.

12.3. Bioaccumulative Potential

Data not available.

12.4. Mobility in Soil

Data not available.

12.5. Other Adverse Ecological Effects

Data not available.

SECTION 13: Disposal Considerations

13.1. Waste Treatment Methods

Data not available.

SECTION 14: Transportation Information

14.1. Transportation by Land-Department of Transportation (DOT, United States of America)

Not regulated according to DOT Regulations.

14.2. Transportation by Air - International Air Transport Association (IATA)

Not regulated according to IATA Dangerous Goods Regulations.

14.3 Transportation of Dangerous Goods (TDG, Canada)

Not regulated according to TDG Regulations.

SECTION 15: Regulatory Information

15.1. Occupational Safety and Health Administration (OSHA) Hazards

Not listed.

15.2. Superfund Amendments and Reauthorization Act (SARA) 302 Extremely Hazardous Substances

Not listed.

15.3. Superfund Amendments and Reauthorization Act (SARA) 311/312 Hazardous Chemicals

Sodium Phosphate Dibasic (CAS # 7558-79-4): 5000 lb final RQ (listed under Sodium phosphate, dibasic); 2270 kg final RQ (listed under Sodium phosphate, dibasic)

Sodium Phosphate Dibasic (CAS # 7558-79-4): 5000 lb final RQ; 2270 kg final RQ

15.4. Superfund Amendments and Reauthorization Act (SARA) 313 Toxic Release Inventory (TRI)

Not listed.

15.5. Massachusetts Right-to-Know Substance List

Sodium Phosphate Dibasic (CAS # 7558-79-4): Present

15.6. Pennsylvania Right-to-Know Hazardous Substances

Sodium Phosphate Dibasic (CAS # 7558-79-4): Environmental hazard

Sodium Phosphate Dibasic (CAS # 7558-79-4): Present

Water (CAS # 7732-18-5): Present

15.7. New Jersey Worker and Community Right-to-Know Components

Sodium Phosphate Dibasic (CAS # 7558-79-4): sn 1723

15.8. California Proposition 65

Not listed.

15.9. Canada Domestic Substances List / Non-Domestic Substances List (DSL/NDSL)

Potassium Chloride (CAS # 7447-40-7): Present (DSL)

Sodium Phosphate Dibasic (CAS # 7558-79-4): Present (DSL)

Sodium Chloride (CAS # 7647-14-5): Present (DSL)

Water (CAS # 7732-18-5): Present (DSL)

Potassium Dihydrogen Phosphate (CAS # 7778-77-0): Present (DSL)

15.10. United States of America Toxic Substances Control Act (TSCA) List

All components of this solution are listed as active on the TSCA Inventory or are mixtures (hydrates) of active items listed on the TSCA Inventory.

Potassium Chloride (CAS # 7447-40-7): Present

Sodium Phosphate Dibasic (CAS # 7558-79-4): Present

Sodium Chloride (CAS # 7647-14-5): Present

Water (CAS # 7732-18-5): Present

Potassium Dihydrogen Phosphate (CAS # 7778-77-0): Present

15.11. European Inventory of Existing Commercial Chemical Substances (EINECS), European

List of Notified Chemical Substances (ELINCS), and No Longer Polymers (NLP)

Potassium Chloride (CAS # 7447-40-7): 231-211-8

Sodium Phosphate Dibasic (CAS # 7558-79-4): 231-448-7

Sodium Phosphate Dibasic (CAS # 7558-79-4): 231-558-5

Sodium Chloride (CAS # 7647-14-5): 231-598-3

Water (CAS # 7732-18-5): 231-791-2

Water (CAS # 7732-18-5): 232-148-9

Potassium Dihydrogen Phosphate (CAS # 7778-77-0): 231-913-4

Potassium Dihydrogen Phosphate (CAS # 7778-77-0): 240-213-8

SECTION 16: Other Information

16.1. Full Text of Hazard Statements and Precautionary Statements

Causes serious eye irritation.

Wash arms, hands and face thoroughly after handling. Wear protective gloves and eye protection.

IF IN EYES: Rinse cautiously with water for several minutes. Remove contact lenses, if present and easy to do. Continue rinsing. If eye irritation persists: Get medical attention.

16.2. Miscellaneous Hazard Classes

Canadian Carcinogenicity Hazard Class: Not Applicable.


Physical Hazards Not Otherwise Classified (PHNOC): Not Applicable.

Health Hazards Not Otherwise Classified (HHNOC): Not Applicable.

Biohazardous Infectious Materials Hazard Class: Not Applicable.

16.3. National Fire Protection Association (NFPA) Rating

Health: 1
Flammability: 0
Reactivity: 0
Special Hazard:



16.4. Document Revision

Last Revision Date: 01/01/2019

DISCLAIMER

When handled properly by qualified personnel, the product described herein does not present a significant health or safety hazard. Alteration of its characteristics by concentration, evaporation, addition of other substances, or other means may present hazards not specifically addressed herein and which must be evaluated by the user. The information furnished herein is believed to be accurate and represents the best data currently available to us. No warranty, expressed or implied, is made and Dawn scientific Inc assumes no legal responsibility or liability whatsoever resulting from its use.