

Product Specification

Product Information:

Products:	Buffer solution pH 1.68
Catalog:	B16013
Grade:	Buffer solution

Reagent Used:

Reagents	CAS #	Grade
Water	7732-18-5	USP/ASTM
Oxalic acid dihydrate	6153-56-6	ACS
Potassium hydroxide	1310-58-3	ACS

Specification:

Test	Specification
Appearance	Clear colorless solution
pH @ 25 °C verified against NIST buffer	pH 1.68 ± 0.01

- For laboratory use only.
- The contents of Product specifications are subject to change.
- Storage: For information refer respective product safety data sheet (SDS).

