

## Product Specification

### Product Information:

<b>Products:</b>	Buffer solution pH 10 at 25°C
<b>Catalog:</b>	B16016
<b>Grade:</b>	Reagent Solution

### Reagent Used:

Reagents	CAS #	Grade
Water	7732-18-5	ASTM/USP
Sodium carbonate	497-19-8	ACS
Sodium bicarbonate	144-55-8	ACS
Sodium hydroxide	1310-73-2	ACS

### Specification:

Test	Specification
Appearance	Clear blue solution
pH @ 25 degree C verified against NIST buffer	pH 10.00 ± 0.01

- For laboratory use only.
- The contents of Product specifications are subject to change.
- Storage: For information refer respective product safety data sheet (SDS).

