

Product Specification

Product Information:

Products:	Buffer solution pH 4.01
Catalog:	B16022
Grade:	Buffer solution

Reagent Used:

Reagents	CAS #	Grade
Water	7732-18-5	USP/ASTM
Formaldehyde	50-00-0	ACS
Potassium hydrogen phthalate	887-24-7	ACS

Specification:

Test	Specification
Appearance	Clear red solution
pH @ 25 °C verified against NIST buffer	pH 4.01 ± 0.01

- For laboratory use only.
- The contents of Product specifications are subject to change.
- Storage: For information refer respective product safety data sheet (SDS).

