

## Product Specification

### Product Information:

<b>Products:</b>	Electrolyte KCl Solution
<b>Catalog:</b>	B16045
<b>Grade:</b>	Buffer solution

### Reagent Used:

Reagents	CAS #	Grade
Water	7732-18-5	USP/ASTM
Potassium chloride	7447-40-7	ACS

### Specification:

Test	Specification
Appearance	Colorless liquid
M of Potassium chloride	3.0 M $\pm$ 0.1M

- For laboratory use only.
- The contents of Product specifications are subject to change.
- Storage: For information refer respective product safety data sheet (SDS).

