

Product Specification

Product Information:

Products:	Sodium Hydroxide 30% W/V
Catalog:	BPS15530
Grade:	Percentage solution

Reagent Used:

Reagents	CAS #	Grade
Water	7732-18-5	USP/ASTM
Sodium Hydroxide	1310-73-2	ACS

Specification:

Test	Specification
Appearance	Colorless Liquid
Assay	30.0 % ± 1.0%

- For laboratory use only.
- The contents of Product specifications are subject to change.
- Storage: For information refer respective product safety data sheet (SDS).

Verified By



Anjal Patel

