

Product Specification

Product Information:

Products:	Sodium Hydroxide 50% w/v
Catalog:	BPS15534
Grade:	Reagent Solution

Reagent Used:

Reagents	CAS #	Grade
Water	7732-18-5	ACS/USP/ASTM
Sodium Hydroxide	1310-73-2	ACS

Specification:

Test	Specification
Appearance	Clear liquid
Assay	49.0 – 51.0 % w/v

- For laboratory use only.
- The contents of Product specifications are subject to change.
- Storage: For information refer respective product safety data sheet (SDS).

Verified By



Anjal Patel

