

Product : Ammonium chloride
Catalogue : CH1031
CAS : 12125-02-9
Formula : NH₄Cl
MW : 53.49
Grade : ACS

Test	Specification
Appearance	White to off white powder/solid.
Assay	≥ 99.5%
Insoluble matter	≤ 0.005%
Sulfate (SO ₄)	≤ 0.002%
Calcium (Ca)	≤ 0.001%
Magnesium (Mg)	≤ 0.0005%
Heavy metals	≤ 0.0005%
Iron (Fe)	≤ 0.0002%
Phosphate (PO ₄)	≤ 0.0002%

- For laboratory use only.
- The contents of Product specifications are subject to change.
- Storage: For information refer respective product safety data sheet (SDS).

