

Product : Ammonium molybdate tetrahydrate
Catalogue : CH1041
CAS : 12054-85-2
Formula : $(\text{NH}_4)_6\text{Mo}_7\text{O}_{24} \cdot 4\text{H}_2\text{O}$
MW : 1236.00
Grade : ACS

Test	Specification
Appearance	Colorless or white solid
Assay	Between 81.0% to 83.0%
Insoluble matter	$\leq 0.005\%$
Nitrate (NO_3)	To pass the test
Chloride (Cl)	$\leq 0.002\%$
Sulfate (SO_4)	$\leq 0.02\%$
Heavy metals	$\leq 0.001\%$
Arsenate, Phosphate and silicate (as SiO_2)	$\leq 0.001\%$
Phosphate (PO_4)	$\leq 0.0005\%$
Magnesium (Mg)	$\leq 0.005\%$
Sodium (Na)	$\leq 0.01\%$
Potassium (K)	$\leq 0.01\%$

- For laboratory use only.
- The contents of Product specifications are subject to change.
- Storage: For information refer respective product safety data sheet (SDS).

