

Product : Ammonium thiocyanate  
Catalogue : CH1055  
CAS : 1762-95-4  
Formula : NH<sub>4</sub>SCN  
MW : 76.12 g/mol  
Grade : -

Test	Specification
Appearance	Colorless or white crystals
Assay	≥ 97.5 %
Residue on ignition	≤ 0.025 %
Insoluble matter	≤ 0.005 %
Chloride Content	≤ 0.005 %
pH of 5 % solution	4.5 to 6.0
Sulfate	≤ 0.005 %
Iron (Fe)	≤ 0.0003 %
Heavy Metals (as Pb)	≤ 0.0005 %

- For laboratory use only.
- The contents of Product specifications are subject to change.
- Storage: For information refer respective product safety data sheet (SDS).

