

Product : Arsenic trioxide
Catalogue : CH1060
CAS : 1327-53-3
Formula : As₂O₃
MW : 197.84
Grade : ACS

Test	Specification
Appearance	White solid.
Assay	99.8 to 100.2%
Residue after ignition	≤ 0.02%
Chloride Content	≤ 0.005 %
Iron (Fe)	≤ 0.0005 %
Insoluble in dilute hydrochloric acid	≤ 0.01% %
Sulfide (S)	To pass the test
Antimony (Sb)	≤ 0.05 %
Lead (Pb)	≤ 0.002 %

- For laboratory use only.
- The contents of Product specifications are subject to change.
- Storage: For information refer respective product safety data sheet (SDS).

