

Product : Barium chloride dihydrate
Catalogue : CH1067
CAS : 10326-27-9
Formula : BaCl₂.2H₂O
MW : 244.26 g/mol
Grade : ACS

Test	Specification
Appearance	White to Off-White crystalline powder/crystals
Assay	≥ 99.0 %
pH @ 25 Deg C (5%, Solution)	5.2 - 8.2
Insoluble matter	≤ 0.005 %
Oxidizing Substances (as NO ₃)	≤ 0.005 %
Loss on Drying	14.0 - 16.0 %
Calcium (Ca)	≤ 0.05 %
Potassium (K)	≤ 0.0025 %
Sodium (Na)	≤ 0.005 %
Strontium (Sr)	≤ 0.1 %
Heavy Metals	≤ 5 ppm
Iron (Fe)	≤ 2 ppm

- For laboratory use only.
- The contents of Product specifications are subject to change.
- Storage: For information refer respective product safety data sheet (SDS).

