

Product : Calcium chloride anhydrous
Catalogue : CH1112
CAS : 10043-52-4
Formula : CaCl₂
MW : 110.98
Grade : ACS/USP

Test	Specification
Appearance	White crystalline powder/solid
Assay	99.0 % to 107.0%
Acid-insoluble matter	≤ 0.02%
Heavy metals	≤ 0.001%
Identification test	To pass the test
Heavy metals	≤ 0.001%
Limit of magnesium and alkali salts	≤ 1.0%
Iron, aluminum and phosphate	To comply the test
Fluoride	≤ 0.004%

- For laboratory use only.
- The contents of Product specifications are subject to change.
- Storage: For information refer respective product safety data sheet (SDS).

