

Product : Calcium hydroxide powder
Catalogue : CH1120
CAS : 1305-62-0
Formula : Ca(OH)₂
MW : 74.09
Grade : ACS

Test	Specification
Appearance	White powder
Assay	≥ 95.0 %
Content of % CaCO ₃	≤ 3.0 %
Insoluble in hydrochloric acid	≤ 0.03 %
Chloride Content	≤ 0.03 %
Iron (Fe)	≤ 0.05 %
Heavy Metals	≤ 0.003 %
Potassium (K)	≤ 0.5 %
Magnesium (Mg)	≤ 0.5 %
Sodium (Na)	≤ 0.05 %
Stronium (Sr)	≤ 0.05%
Sulfur compounds (as SO ₄)	≤ 0.10 %

- For laboratory use only.
- The contents of Product specifications are subject to change.
- Storage: For information refer respective product safety data sheet (SDS).

