

Product : Cupric sulfate pentahydrate
Catalogue : CH1195
CAS : 7758-99-8
Formula : $\text{CuSO}_4 \cdot 5\text{H}_2\text{O}$
MW : 249.69 g/mol
Grade : ACS

Test	Specification
Appearance	Blue Crystals/solid
Assay	98.0 - 102.0 %
Insoluble matter	≤ 0.005 %
Chloride Content	≤ 0.001 %
Iron (Fe)	≤ 0.003 %
Nitrogen Compounds	≤ 0.002 %
Calcium (Ca)	≤ 0.005 %
Potassium (K)	≤ 0.01 %
Sodium (Na)	≤ 0.02 %
Nickel (Ni)	≤ 0.005 %

- For laboratory use only.
- The contents of Product specifications are subject to change.
- Storage: For information refer respective product safety data sheet (SDS).

