

Product : EDTA disodium salt dihydrate
Catalogue : CH1209
CAS : 6381-92-6
Formula : $C_{10}H_{14}N_2O_8Na_2 \cdot 2H_2O$
MW : 372.24 g/mol
Grade : ACS

Test	Specification
Appearance	White powder or crystals or solid.
Assay	99.0% to 101.0%
pH of 5% solution at 25°C	4.0 to 6.0
Insoluble matter	≤ 0.005%
Heavy metals (As Pb)	≤ 0.005%
Nitilotriacetic acid	≤ 0.1 %
Iron (Fe)	≤ 0.01 %

- For laboratory use only.
- The contents of Product specifications are subject to change.
- Storage: For information refer respective product safety data sheet (SDS).

