

Product : Potassium hydrogen phthalate
Catalogue : CH1382
CAS : 877-24-7
Formula : HOCC6H4COOK
MW : 204.22 g/mol
Grade : ACS

Test	Specification
Appearance	White Crystals
Assay (on dried basis)	99.5 - 100.5 %
Insoluble matter	≤ 0.005 %
Chlorine compounds (Cl)	≤ 0.003 %
Iron (Fe)	≤ 5 ppm
Heavy metal (as Pb)	≤ 5 ppm
Sodium (Na)	≤ 0.005 %
Sulfur compounds (as S)	≤ 0.002 %

- For laboratory use only.
- The contents of Product specifications are subject to change.
- Storage: For information refer respective product safety data sheet (SDS).

