

Product : Sodium acetate anhydrous
Catalogue : CH1494
CAS : 127-09-3
Formula : $C_2H_3NaO_2$
MW : 82.03 g/mol
Grade : ACS

Test	Specification
Appearance	White powder or crystals
Assay	≥ 99.0 %
Loss on Drying	≤ 1.0 %
Insoluble matter	≤ 0.01 %
Chloride Content	≤ 0.002 %
pH of 5 % solution	7.0 - 9.2
Phosphate (PO ₄)	≤ 0.001 %
Sulfate	≤ 0.003 %
Iron (Fe)	≤ 0.001 %
Heavy Metals	≤ 0.001 %
Calcium (Ca)	≤ 0.005 %
Magnesium (Mg)	≤ 0.002 %

- For laboratory use only.
- The contents of Product specifications are subject to change.
- Storage: For information refer respective product safety data sheet (SDS).

