

Product : Potassium chlorate  
Catalogue : CH1504  
CAS : 3811-04-9  
Formula :  $\text{KClO}_3$   
MW : 122.55 g/mol  
Grade : Lab

Test	Specification
Appearance	White crystalline/granular solid.
Assay	$\geq 98.0\%$
Insoluble matter	$\leq 0.005\%$
Bromate	$\leq 0.015\%$
Chloride content	$\leq 0.001\%$
Sulfate	To pass the test
Iron (Fe)	$\leq 3\text{ppm}$
Magnesium (Mg)	$\leq 0.002\%$
Sodium (Na)	$\leq 0.01\%$

- For laboratory use only.
- The contents of Product specifications are subject to change.
- Storage: For information refer respective product safety data sheet (SDS).

