

Product : Sodium carbonate anhydrous
Catalogue : CH1520
CAS : 497-19-8
Formula : Na_2CO_3
MW : 105.89
Grade : ACS

Test	Specification
Appearance	Colorless to white solid/Powder.
Assay	$\geq 99.5\%$
Loss on heating at 285°C	$\leq 1.0\%$
Silicate	$\leq 0.005\%$
Chloride (Cl)	$\leq 0.001\%$
Sulfate (SO_4)	$\leq 0.003\%$
Heavy metals (as Pb)	$\leq 0.001\%$
Iron (Fe)	$\leq 0.001\%$
Phosphate (PO_4)	$\leq 0.001\%$
Calcium (Ca)	$\leq 0.03\%$
Magnesium (Mg)	$\leq 0.005\%$
Potassium (K)	$\leq 0.005\%$

- For laboratory use only.
- The contents of Product specifications are subject to change.
- Storage: For information refer respective product safety data sheet (SDS).

