

Product : Sodium citrate dihydrate
Catalogue : CH1536
CAS : 6132-04-3
Formula : $C_6H_5Na_3O_7 \cdot 2H_2O$
MW : 294.10 g/mol
Grade : ACS

Test	Specification
Appearance	White powder or crystals
Odor	Odorless
Assay	≥ 99.0%
pH of 5% solution	7.00 to 9.00
Alkalinity	To pass the test
Insoluble matter	≤ 0.005%
Heavy metals	≤ 5 ppm
Chloride (Cl)	≤ 0.003 %
Sulfate (SO ₄)	≤ 0.005 %
Calcium (Ca)	≤ 0.005 %
Ammonia (NH ₃)	≤ 0.003 %
Iron (Fe)	≤ 5 ppm

- For laboratory use only.
- The contents of Product specifications are subject to change.
- Storage: For information refer respective product safety data sheet (SDS).

