

Product : Sodium hydroxide pellet  
Catalogue : CH1560  
CAS : 1310-73-2  
Formula : NaOH  
MW : 40.00 g/mol  
Grade : ACS

Test	Specification
Appearance	White Pellets
Assay	≥ 97.0 %
Sodium carbonate	≤ 1.0 %
Sulfate (SO <sub>4</sub> )	≤ 0.003 %
Chloride Content	≤ 0.005 %
Nitrogen Compounds (as N)	≤ 0.001 %
Phosphate (PO <sub>4</sub> )	≤ 0.001 %
Heavy Metals (as Ag)	≤ 0.002 %
Iron (Fe)	≤ 0.001 %
Nickel (Ni)	≤ 0.001 %
Mercury (Hg)	≤ 0.1 ppm
Calcium (Ca)	≤ 0.005 %
Magnesium (Mg)	≤ 0.002 %
Potassium (K)	≤ 0.02 %

- For laboratory use only.
- The contents of Product specifications are subject to change.
- Storage: For information refer respective product safety data sheet (SDS).

