

Product : Sodium Phosphate dibasic anhydrous
Catalogue : CH1589
CAS : 7558-79-4
Formula : Na_2HPO_4
MW : 141.96 g/mol
Grade : ACS

Test	Specification
Appearance	White to Off White crystalline powder.
Assay	$\geq 99.0\%$
pH 5% solution at 25°C	8.7 to 9.3
Insoluble Matter	$\leq 0.01\%$
Loss on Drying	$\leq 0.2\%$
Chloride (Cl)	$\leq 0.002\%$
Sulfate (SO ₄)	$\leq 0.005\%$
Heavy Metals	$\leq 0.001\%$
Iron (Fe)	$\leq 0.002\%$

- For laboratory use only.
- The contents of Product specifications are subject to change.
- Storage: For information refer respective product safety data sheet (SDS).

