

Product : Sodium phosphate tribasic dodecahydrate
Catalogue : CH1595
CAS : 10101-89-0
Formula : $\text{Na}_3\text{PO}_4 \cdot 12\text{H}_2\text{O}$
MW : 380.12 g/mol
Grade : ACS

Test	Specification
Appearance	Colorless to white solid.
Assay	Between 98.0 to 102.0%
Excess alkali (as NaOH)	$\leq 2.5\%$
Insoluble matter	$\leq 0.01\%$
Chloride (Cl)	$\leq 0.001\%$
Sulfate (SO_4)	$\leq 0.01\%$
Heavy metals (as Pb)	$\leq 0.001\%$
Iron (Fe)	$\leq 0.001\%$

- For laboratory use only.
- The contents of Product specifications are subject to change.
- Storage: For information refer respective product safety data sheet (SDS).

