

Product : Sodium sulfate anhydrous
Catalogue : CH1604
CAS : 7757-82-6
Formula : Na₂SO₄
MW : 142.04 g/mol
Grade : ACS

Test	Specification
Appearance	White crystalline/granular solid.
Assay	≥ 99.0%
pH (5 % solution at 25°C)	5.2 to 9.2
Insoluble matter	≤ 0.001 %
Loss on ignition	≤ 0.5%
Chloride (Cl)	≤ 0.001 %
Nitrogen compounds	≤ 5ppm
Calcium (Ca)	≤ 0.01 %
Heavy metals (by ICP-OES)	≤ 5ppm
Iron (Fe)	≤ 0.001%
Phosphate (PO ₄)	≤ 0.001%
Potassium (K)	≤ 0.01%
Magnesium (Mg)	≤ 0.005%

- For laboratory use only.
- The contents of Product specifications are subject to change.
- Storage: For information refer respective product safety data sheet (SDS).

