

Product : Titanium dioxide
Catalogue : CH1659
CAS : 13463-67-7
Formula : TiO₂
MW : 79.87
Grade : ACS

Test	Specification
Appearance	White to almost white powder to crystalline
Purity	Min. 98.0%
Residue on ignition (1000°C for 2 Hr)	≤ 0.5%
Arsenic (As)	≤ 0.002%
Iron (Fe)	≤ 0.005%
Lead (Pb)	≤ 0.005%

- For laboratory use only.
- The contents of Product specifications are subject to change.
- Storage: For information refer respective product safety data sheet (SDS).

