

Product : Acetone
Catalogue : CH5001
CAS : 67-64-1
Formula : CH₃COCH₃
MW : 58.08
Grade : ACS

Test	Specification
Appearance	Colorless liquid.
Color (IN APHA)	≤ 10 APHA
Density (g/ml) at 20°C	0.791 to 0.793
Titrateable acid	≤ 0.0003 meq/g
Titrateable base	≤ 0.0006 meq/g
Assay by GC	≥ 99.5 %
Methanol	≤ 0.05 %
Isopropyl alcohol	≤ 0.05 %
Solubility	To pass the test
Water	≤ 0.5 %
Residue on evaporation	≤ 0.001 %
Substance reducing substances	To pass the test
Aldehydes (as HCHO)	≤ 0.002 %

- For laboratory use only.
- The contents of Product specifications are subject to change.
- Storage: For information refer respective product safety data sheet (SDS).

