

Product : Acetone Lab
Catalog# : CH5002
CAS : 67-64-1
Formula : CH₃COCH₃
MW : 58.08

TEST	SPECIFICATION
Appearance	Colorless liquid.
Color (IN APHA)	≤ 10 APHA
Titrateable acid	≤ 0.003 meq/g
Titrateable base	≤ 0.006 meq/g
Assay by GC	≥ 99.0 %
Methanol	≤ 0.05 %
Isopropyl alcohol	≤ 0.05 %
Solubility	To pass the test
Water	≤ 0.5 %
Residue on evaporation	≤ 0.01 %

- For laboratory use only.
- The contents of Product specifications are subject to change.
- Storage: For information refer respective product safety data sheet (SDS).

