

Product : Acetonitrile ACS
Catalog# : CH5003
CAS : 75-05-8
Formula : CH₃CN
MW : 41.05
Grade : ACS

TEST	SPECIFICATION
Appearance	Clear, colorless liquid
Assay	≥ 99.5%
Color	≤ 10 APHA
Residue after evaporation	≤ 0.005%
Titration acid	≤ 8 meq/g
Titration base	≤ 0.6 meq/g
Water	≤ 0.3%

- For laboratory use only.
- The contents of Product specifications are subject to change.
- Storage: For information refer respective product safety data sheet (SDS).

