

Product : Benzene ACS
Catalog# : CH5005
CAS : 71-43-2
Formula : C₆H₆
MW : 78.11
Grade : ACS

TEST	SPECIFICATION
Appearance	Clear liquid
Assay	≥ 99.0 %
Color (APHA)	≤ 10
Residue after evaporation	≤ 0.001 %
Substance darkened by sulfuric acid	Passes test
Thiophene (limit about 1 ppm)	Passes test
Sulfur compound (as S)	≤ 0.005 %
Water (H ₂ O)	≤ 0.05 %

- For laboratory use only.
- The contents of Product specifications are subject to change.
- Storage: For information refer respective product safety data sheet (SDS).

