

Product : Benzyl Alcohol ACS
Catalog# : CH5007
CAS : 100-51-6
Formula : C₆H₅CH₂OH
MW : 108.14
Grade : ACS

TEST	SPECIFICATION
Appearance	Clear, colorless liquid
Assay	≥ 99.0 %
Color (APHA)	≤ 20
Residue after ignition	≤ 0.005%
Acetophenone (C ₆ H ₅ COCH ₃)	≤ 0.02 %
Benzaldehyde (C ₆ H ₅ CHO)	≤ 0.01%

- For laboratory use only.
- The contents of Product specifications are subject to change.
- Storage: For information refer respective product safety data sheet (SDS).

