

Product : Benzyl alcohol FCC  
Catalog# : CH5008  
CAS : 100-51-6  
Formula : C<sub>7</sub>H<sub>8</sub>O  
MW : 108.14  
Grade : FCC

TEST	SPECIFICATION
Appearance	Colorless liquid with sharp burning taste and faint aromatic odor
Assay	≥ 99.0 %
Identification	Conforms to Infrared Spectra
Refractive Index	1.539 – 1.541 @ 20C
Specific Gravity	1.042 – 1.047 @ 25C
Distillation Range	≤ 95% distills between 202.5-206.5C
Aldehydes	≤ 0.2%
Chlorinated Compounds	To Pass Test

- For laboratory use only.
- The contents of Product specifications are subject to change.
- Storage: For information refer respective product safety data sheet (SDS).

