

Product : n-Butyl alcohol
Catalogue : CH5010
CAS : 71-36-3
Formula : $\text{CH}_3(\text{CH}_2)_2\text{CH}_2\text{OH}$
MW : 74.12 g/mol
Grade : ACS

Test	Specification
Appearance	Colorless liquid.
Color (APHA)	≤ 10 APHA
Purity (GC)	$\geq 99.40\%$
Butyl ether	$\leq 0.2\%$
Carbonyl compounds (as Butyraldehyde)	$\leq 0.01\%$
Residue on evaporation	$\leq 0.005\%$
Water (by Karl fisher)	$\leq 0.1\%$
Titration acid (meq/g)	≤ 0.0008 meq/g

- For laboratory use only.
- The contents of Product specifications are subject to change.
- Storage: For information refer respective product safety data sheet (SDS).

