

Product : tert-Butyl alcohol ACS
Catalog# : CH5011
CAS : 75-65-0
Formula : C₄H₁₀O
MW : 74.12
Grade : ACS

TEST	SPECIFICATION
Appearance	Colorless liquid or solid.
Purity (GC)	≥ 99.0%
Identification	Conforms to Infrared Spectra
Water	≤ 0.1%
Residue on evaporation	≤ 0.003%
Color test (APHA)	≤ 20 APHA
Carbonyl compounds as formaldehyde	≤ 0.01%
Titrateable acid (meq/g)	≤ 0.001 meq/g

- For laboratory use only.
- The contents of Product specifications are subject to change.
- Storage: For information refer respective product safety data sheet (SDS).

