

Product : Chlorobenzene ACS
Catalog# : CH5012
CAS : 108-90-7
Formula : C₆H₅Cl
MW : 112.56
Grade : ACS

TEST	SPECIFICATION
Appearance	Clear liquid
Assay	≥ 99.5 %
Color	≤ 30
Residue after evaporation	≤ 0.02 %
Titration acid	≤ 0.004 meq/g

- For laboratory use only.
- The contents of Product specifications are subject to change.
- Storage: For information refer respective product safety data sheet (SDS).

