

Product : Chloroform
Catalogue : CH5013
CAS : 67-66-3
Formula : CHCl_3
MW : 119.38 g/mol
Grade : ACS

Test	Specification
Appearance	Clear colorless liquid.
Color	≤ 10 APHA
Acetone and aldehyde	To pass the test
Purity (GC)	$\geq 99.8\%$
Acid and chloride	To pass the test
Lead (Pb)	≤ 0.05 ppm
Substances darkened by sulfuric acid	To pass the test
Residue after evaporation	$\leq 0.001\%$
Free chloride (Cl)	To pass the test
Water (By KF)	$\leq 0.03\%$

- For laboratory use only.
- The contents of Product specifications are subject to change.
- Storage: For information refer respective product safety data sheet (SDS).

