

Product : Cyclohexane  
Catalogue : CH5016  
CAS : 110-82-7  
Formula : C<sub>6</sub>H<sub>12</sub>  
MW : 84.16  
Grade : ACS

Test	Specification
Appearance	Clear colorless liquid
Color Test (APHA)	≤ 10 APHA
Purity (GC)	≥ 99.0 %
Residue on Evaporation	≤ 0.002 %
Water (by Karl Fischer)	≤ 0.02%
Substances Darkened (by H <sub>2</sub> SO <sub>4</sub> )	To pass the test

- For laboratory use only.
- The contents of Product specifications are subject to change.
- Storage: For information refer respective product safety data sheet (SDS).

