

Product : Cyclohexane Reagent  
Catalog# : CH5017  
CAS : 110-82-7  
Formula : C<sub>6</sub>H<sub>12</sub>  
MW : 84.16

TEST	SPECIFICATION
Appearance	Clear colorless liquid
Color Test (APHA)	≤ 10 APHA
Purity (GC)	≥ 98.5 %
Residue on Evaporation	≤ 0.002 %
Water (by Karl Fischer)	≤ 0.02%

- For laboratory use only.
- The contents of Product specifications are subject to change.
- Storage: For information refer respective product safety data sheet (SDS).

