

Product : Cyclohexanone ACS  
Catalog# : CH5018  
CAS : 108-94-1  
Formula : C<sub>6</sub>H<sub>10</sub>O  
MW : 98.14  
Grade : ACS

TEST	SPECIFICATION
Appearance	Clear, oily liquid
Assay	≥ 99.0 %
Color (APHA)	≤ 10
Residue after evaporation	≤ 0.05 %
Water (H <sub>2</sub> O)	≤ 0.05 %

- For laboratory use only.
- The contents of Product specifications are subject to change.
- Storage: For information refer respective product safety data sheet (SDS).

