

Product : Decane ACS
Catalog# : CH5020
CAS : 124-18-5
Formula : $\text{CH}_3(\text{CH}_2)_8\text{CH}_3$
MW : 142.28
Grade : ACS

TEST	SPECIFICATION
Appearance (Colour)	Colorless Liquid
Infrared Spectrum	Conforms to Structure
Refractive index at 20°C	1.409 - 1.413
Purity (GC)	≥ 99.0 %

- For laboratory use only.
- The contents of Product specifications are subject to change.
- Storage: For information refer respective product safety data sheet (SDS).

