

Product : 1,2-Dichloroethane ACS
Catalog# : CH5022
CAS : 107-06-2
Formula : CH₂ClCH₂Cl
MW : 98.95
Grade : ACS

TEST	SPECIFICATION
Appearance	Clear, colorless liquid
Assay	≥ 99.0 %
Color (APHA)	≤ 10
Residue after evaporation	≤ 0.002 %
Titration acid	≤ 0.0003 %
Water (H ₂ O)	≤ 0.03 %

- For laboratory use only.
- The contents of Product specifications are subject to change.
- Storage: For information refer respective product safety data sheet (SDS).

