

Product : Dichloromethane (Methylene chloride)
Catalogue : CH5024
CAS : 75-09-2
Formula : CH₂Cl₂
MW : 84.93 g/mol
Grade : ACS

Test	Specification
Appearance	Colorless Liquid
Color (APHA)	≤ 10 APHA
Purity (GC area %)	≥ 99.80 %
Density (g/ml) at 25°C	1.315 to 1.321
Water (By KF)	≤ 0.02%
Residue on evaporation	≤ 0.002 %
Titrateable acid (meq/g)	≤ 0.0003 meq/g
Free halogens	To pass the test

- For laboratory use only.
- The contents of Product specifications are subject to change.
- Storage: For information refer respective product safety data sheet (SDS).

