

Product : Dichloromethane (Methylene Chloride) Lab
Catalog# : CH5025
CAS : 75-09-2
Formula : CH₂Cl₂
MW : 84.93 g/mol

TEST	SPECIFICATION
Appearance	Colorless Liquid
Color (APHA)	≤ 10 APHA
Purity (GC area %)	≥ 99.80 %
Water (By KF)	≤ 0.02%
Residue on evaporation	≤ 0.02 %
Titrateable acid (meq/g)	≤ 0.003 meq/g

- For laboratory use only.
- The contents of Product specifications are subject to change.
- Storage: For information refer respective product safety data sheet (SDS).

