

Product : Diethyl Ether
Catalogue : CH5026
CAS : 60-29-7
Formula : (CH₃CH₂)₂O
MW : 74.12
Grade : ACS

Test	Specification
Appearance (color)	Colorless
Appearance (form)	Liquid
Purity (GC area %)	≥ 99.0 %
Water	≤ 0.03%
Color (APHA)	≤ 10 APHA
Residue after evaporation	≤ 0.001%
Titration acid	≤ 0.0002 meq/g
Carbonyl (As HCHO)	≤ 0.001%
Alcohol (CH ₂ CH ₂ OH)	To pass the test
Stabilizer (BHT)	Between 5 to 10 ppm

- For laboratory use only.
- The contents of Product specifications are subject to change.
- Storage: For information refer respective product safety data sheet (SDS).

