

Product : Diethyl Ether Anhydrous ACS
Catalog# : CH5027
CAS : 60-29-7
Formula : C₄H₁₀O
MW : 74.12
Grade : ACS

TEST	SPECIFICATION
Appearance	Clear colorless liquid
Assay	≥ 99.0 %
Color (APHA)	≤ 10
Peroxide (as H ₂ O ₂)	≤ 1 ppm
Residue after evaporation	≤ 0.001 %
Titration acid	≤ 0.0002 meq/g
Carbonyl (as HCHO)	≤ 0.001 %
Alcohol (CH ₃ CH ₂ OH)	To pass the test
Water (H ₂ O)	≤ 0.03 %

- For laboratory use only.
- The contents of Product specifications are subject to change.
- Storage: For information refer respective product safety data sheet (SDS).

