

Product : Dimethyl Sulfoxide ACS
Catalog# : CH5029
CAS : 67-68-5
Formula : $(\text{CH}_3)_2\text{SO}$
MW : 78.13
Grade : ACS

TEST	SPECIFICATION
Appearance	Clear, colorless liquid
Assay	$\geq 99.9\%$
Residue after evaporation	0.01 %
Titration acid	0.001 meq/g
Water (H ₂ O)	0.01 %

- For laboratory use only.
- The contents of Product specifications are subject to change.
- Storage: For information refer respective product safety data sheet (SDS).

