

Product : Dimethylformamide Reagent
Catalog# : CH5031
CAS : 68-12-2
Formula : $(\text{CH}_3)_2\text{NC}(\text{O})\text{H}$
MW : 99.8

TEST	SPECIFICATION
Assay	$\geq 99.8\%$ min.
Color (APHA)	≤ 10
Residue After Evaporation	$\leq 0.005\%$
Titration Base	≤ 0.003 meq/g
Titration Acid	≤ 0.0005 meq/g
Water	$\leq 0.15\%$

- For laboratory use only.
- The contents of Product specifications are subject to change.
- Storage: For information refer respective product safety data sheet (SDS).

