

Product : 1,4 Dioxane ACS  
Catalog# : CH5033  
CAS : 123-91-1  
Formula : C<sub>4</sub>H<sub>8</sub>O<sub>2</sub>  
MW : 88.11  
Grade : ACS

TEST	SPECIFICATION
Appearance	Clear liquid
Assa/gy	≥ 99.0 %
Freezing point	≤ 11.0°C
Color (APHA)	≤ 20
Peroxide (as H <sub>2</sub> O <sub>2</sub> )	≤ 0.005 %
Residue after evaporation	≤ 0.005 %
Titration acid	≤ 0.0016 meq/g
Carbonyl (as HCHO)	≤ 0.01 %
Water (H <sub>2</sub> O)	≤ 0.05 %

- For laboratory use only.
- The contents of Product specifications are subject to change.
- Storage: For information refer respective product safety data sheet (SDS).

