

Product : Ethyl acetate  
Catalogue : CH5034  
CAS : 141-78-6  
Formula :  $\text{CH}_3\text{COOC}_2\text{H}_5$   
MW : 88.11  
Grade : ACS

Test	Specification
Appearance	Colorless Liquid
Purity (GC)	$\geq 99.50\%$
Water (by Karl Fischer)	$\leq 0.2\%$
Residue on Evaporation	$\leq 0.003\%$
Color Test	$\leq 10$ APHA
Substances Darkened (by $\text{H}_2\text{SO}_4$ )	To Pass the test
Titrateable Acid (meq/g)	$\leq 0.0009$
Boiling range	76.2°C to 78.2°C

- For laboratory use only.
- The contents of Product specifications are subject to change.
- Storage: For information refer respective product safety data sheet (SDS).

